



ISOKINETIK METHOD 31

Method 31 SuperVOST Sampling



The Apex Instruments SuperVOST Sampling Kit includes the multifunctional VersaCase2 cabinet with a redesigned mounting panel, glass-lined probe assembly, VOST glassware with clamps, insulated coolant reservoir with a submersible pump and VU-cord. The VU-cord can be used for connection directly to the Model 623 Source Sampler Console, or to Model 572 with the VUA orifice flow adapter.



Call for details

+62 821 4486 1199

For More Information

sales02@totalenviro.co.id

Visit Us

totalenviro.co.id