

ISOKINETIK SYSTEM



Isokinetic sampling system provides every essential piece of apparatus needed to perform Method 5, Method 17, Method 23 runs with one console and additional part of method.

All of the ESC isokinetic sampling systems consist of the following primary components: nozzle set, heated probe assembly with type "S" pitot and stack thermocouple, modular sample case with heated filter box and chilled impinger box, filter and glassware set, impinger outlet adapter with thermocouple, umbilical cable, metering console, rotary vane vacuum pump and suspension rail system. The more advanced systems include a calibration orifice set, transport cases, professional multiple glassware sets, a cleanup and recovery kit, and probes and umbilical cables in a variety of lengths and styles



Call for details

+62 821 4486 1199

For More Information

sales02@totalenviro.co.id

Visit Us

totalenviro.co.id